

Figure 1. A large, ornate, multi-tiered chandelier hanging from a ceiling.

Figure 2. A large, ornate, multi-tiered chandelier hanging from a ceiling.

Figure 3. A large, ornate, multi-tiered chandelier hanging from a ceiling.

Figure 4. A large, ornate, multi-tiered chandelier hanging from a ceiling.

Figure 5. A large, ornate, multi-tiered chandelier hanging from a ceiling.

Figure 6. A large, ornate, multi-tiered chandelier hanging from a ceiling.

Technical Specification Datasheet & Hardware Architecture Overview

The Maber Research Robot is a compelling addition to the world of Research and Education. With the same industry-based hardware of Rethink Robotics' flagship robot Maber that is revolutionizing the manufacturing world, the Maber Research Robot is providing a safe, affordable, robust platform for universities and labs to utilize in education. With market-leading value, a unique inherent safety system and the real-world relevance of the Maber manufacturing solution, it is quickly becoming a must-have tool for leading institutions around the globe.

Maber consists of two, 7-degree-of-freedom arms with five electric actuators at each joint. In comparison, full position and force sensing. Three integrated cameras, along with force, accelerometers and range-finding sensors, Maber's Inuitive Zero Force Gravity Compensation mode allows users to effortlessly move each 7-degree-of-freedom arm and teach the robot positions and trajectories without the need to program them mathematically, enabling collaborative human-robot co-work and advanced Human Robot Interaction.

The Maber Research Robot allows direct programming access to the system via a standard, open-source ROS API interface. Users can run custom programs from a connected development workstation, or locally through access to the on-board CPU. It is entirely safe to operate around humans without the need for safety cages or other guarding equipment, and comes with a suite of example programs, demonstration videos, online documentation and a robust user community to help new users get started quickly in developing new applications.

Further technical details can be found at: <http://dx.doi.org/10.1016/j.robot.2014.06.001>

Open-Source ROS Code can be found at: <https://github.com/RethinkRobotics>

API Documentation can be found at: <http://api.rtrb.com/>

Hardware Architecture Overview

Arms

Robot has two seven-degree-of-freedom (DOF) arms. Seven DOF arms are desirable as they provide a kinematic redundancy greatly improving manipulation flexibility. Each joint contains a Series Elastic Actuator (SEA), this actuation technique is key to making Robot safe.

Series Elastic Actuators differ from traditional actuation techniques in that we introduce a spring between the motor/actuator elements and the output of the actuator. The result is:

- Spring in parallel, compliance force control
- protection against physical loads

Inherent safety is a characteristic of SEA's. The springs in these actuators are deformable by human level inputs. This deformation is an inherent safety mechanism.

By introducing the spring element we are now able to measure torque output from the actuator. This enables such things as torque control, impedance control, and more.

Reinhold Robotic, MIT with the Harvard Paper on the SEA and Robot safety by and at: https://www.csail.mit.edu/pubs/03SEA_robot_safety.html



Robot Movement can use Denavit-Hartenberg's joints and is based using the following convention,
From the Base, along the arm to the End-Effector:

- 00 - Base Roll
- 01 - Base Pitch
- 02 - Elbow Roll
- 03 - Elbow Pitch
- 04 - Wrist Roll
- 05 - Wrist Pitch
- 06 - Wrist Roll



Robot Control

Fundamental to the operation of the robot is the SOE is a hierarchy with the level of joint control modes. The joint control modes define the 'mode' in which the user can control the robot's motion in joint space. The user are provided with a PC software app.

When a joint position command is published from the development PC, the control software processes and represents a robot control plugin according to the message. This message is then parsed, and represented in memory based on control mode.


Depending on the control mode, modifications are made to the input command. These modifications are typically due to the safety control layer (e.g. arm-to-arm collision avoidance, collision detection, etc.)

A control rate timeout is also enforced at the robot control layer. This allows that if a new 'control command' message is not received within the specified timeout (0.3 seconds, or 300ms), the robot will 'Timeout'. When the robot 'Times out', the current control mode is exited, reverting back to position control mode where the robot will command 0rad/s to current joint angles. The reason for this is safety. For example, if controlling in velocity control mode where you are commanding 3.0 Rad/s to joint 01 and you lose network connectivity by the robot could reach a dangerous position. By 'timing out', the robot will be safer, reacting to network timeouts, or incorrect control behaviour.

Joint Position Control

Joint position control mode is the fundamental, basic control mode for Member arm motion. In position control mode, we specify joint angles at which we want the joints to achieve. Typically this will consist of seven values, a commanded position for each of the seven joints, resulting in a full description of the arm configuration. `JointCommand` message: `space_msgs/JointCommand`

This joint command is then subscribed to by Member's Motor Controller. The motor controller processes this joint command and will regulate (collaboration avoidance and default on) and executed behavior by making the following control actions:




Free Joint Position Control

Free joint position control mode was introduced to provide a much more direct position control method. The same best cases as the standard joint position control mode still hold. However, the joint commands are largely left to the JCs for execution.

JointCommandResponse was replaced with

Advanced Control Mode. It no longer apply collision avoidance offsets, and allow for very fast reactions at the joint velocity limit. This mode should be used with care.

The free joint position commands only performs the following special call on typically these solutions are normal:



JointCommandResponse
- JointCommandResponse
- JointCommandResponse
- JointCommandResponse
- JointCommandResponse

Advanced Control Mode
- Advanced Control Mode
- Advanced Control Mode
- Advanced Control Mode

JointCommandResponse
- JointCommandResponse
- JointCommandResponse
- JointCommandResponse

Joint Velocity of the Control


Joint velocity control mode is an advanced control mode. In velocity control mode, we specify joint velocities at which we want the joints to simultaneously achieve. Typical trajectory will be consist of several values, a commanded velocity for each of the several joints.

JointCommandTypeMode: `VELOCITY_CONTROL_MODE`

When commanding joint velocity, if a commanded velocity to one of the joints will result in a joint position that is beyond the joint limits, no joints will be commanded, as all of the commands is considered invalid. Please note that even example of a common control method, using the Jacobian for Cartesian control resulting in joint velocity commands, if a single joint hits its limit, the rest of the joints will still be commanded resulting in obscure and potentially dangerous motions. We recommend either inputting a joint space potential field or joint limit check before submitting the joint velocity commands.

Advanced Control Mode: As well as for very fast joint velocity up to the joint velocity limit, this mode should be used with care.

This joint commands is subscribed to by Member's Motor Controller. The motor controller processes this joint velocity command (applying collision avoidance and detect and expected behavior) by sending the following request:



Joint Torque Control

Joint torque control mode is an advanced control mode. In torque control mode, we specify joint torques at which we want the joints to stand horizontally each time. Typically this will be a set of seven values, a commanded torque for each of the seven joints.


Joint commands during no-op: `OPPOSITE_MPCOM`

Joint torque commands are applied in addition to gravity and spring compensation torques. This default can be disabled by publishing an `std_msgs/String` message on the topic:

```
/robot/1/acknowledge/suppress_gravity_compensation: set  
value = 0.0
```

Advanced Control Mode. When commanding in torque control mode, access is granted to the lowest control levels. This puts much responsibility on the control program and should be used with care.

The joint command is then as described by `Member Motor Controller`. The motor controller processes this joint torque command (applying coil on avoidance and detent on) and expects to stand or by reading the Hall sensor read output.



Motor Torque Command
The motor controller will
stand on the motor
and the motor will
stand on the motor

Motor Position Command
The motor controller will
stand on the motor
and the motor will
stand on the motor

Gravity Compensation Torques

In order to oppose the effect of gravitational force acting on one link of a arm, gravity compensation torques have to be applied across that arm. This is a local torque that is active from the time the robot is enabled. The built-in gravity compensation model uses `JDL` for calculating the gravity compensation torques. The springs attached to `DL` also require compensation torques. These are found using an internal spring model and applied in conjunction with the torques from `JDL`. These gravity compensation torques are passed directly to the `JDL` through a separate channel and are applied in opposite of the control torques.

```
The gravity compensation torques are available via the topic:  
robot/1.1.0.0/robot_controller/parameters/0.0.0.0/0.0.0.0  
message type: builtin_interfaces/msg/Float64Array
```

In order to disable the gravity compensation torques, an empty message should be published to the topic:

```
robot/1.1.0.0/robot_controller/parameters/0.0.0.0/0.0.0.0  
greater than 5 Hz.
```

Zero-G

Zero-G mode can often be confused with the no-deceleration by disabling the gravity compensation torques. By default, the gravity compensation torques will always be applied when the robot is enabled. In Zero-G mode, the controllers are disabled and so the arm can be freely moved across the workspace. In this case, the effect of gravity would be compensated by the gravity compensation model applying a gravity compensation torque across the joints. There would be no torque from the controllers since they would not be active, and so the arm can be moved freely around, hence the name. The Zero-G mode can be enabled by pressing the `0` key on the remote control.



Collision Model

There are a total of three levels of collision models in Rethink.

Each of Rethink's links has its own collision block which is slightly bigger than the link of that link. When these blocks come into contact with each other, collision avoidance model is triggered.

The collision blocks can be visualized in RVIZ by setting the Robot Description Field under RobotModel as

```
~/catkin_ws/src/rethink_robot_description/launch/robot.launch? and enabling the Field Collision Enabled.
```

The first collision model involves the detection of two types of collisions - Impact and Squish.

Impact is detected when a sudden change in torque is sensed across any of the joints. This can be related to a scenario in which a robot might collide with a human. Here, the sudden change in torque is sensed during the impact and the robot comes to a stop for 3 seconds before attempting to move again.

Squish is detected when a joint tries to press against a stationary object. For instance, when a robot arm tries to push a wall, the torque applied across the joint increases and it immediately stops when the applied torque is greater than a threshold. It resumes its motion after 3 seconds.

The third collision model is high speed scaling. This is to avoid the collisions between two links moving at higher speeds.


The collision blocks for this purpose are larger than the original collision blocks.

The collision detection happens when the speed of the link in Cartesian space is higher than 0.3 m/s.

In order to avoid the collisions, the corresponding velocities are scaled down and converted back to position commands.


To disable the collision avoidance on a given arm you need to publish an empty message at greater than 5Hz to the topic

```
~/catkin_ws/src/rethink_robot_description/launch/robot.launch? and enabling the Field Collision Enabled.
```




Software Architecture Overview:

The Member Research Robot (MRR) provides a software interface allowing researchers of all disciplines to develop custom applications to run on the MRR platform. The MRR interface with the Member Research Robot via [ROS \(Robot Operating System\)](#). MRR provides a stand-alone ROS Master to which any development workstation can connect and control MRR via the various [ROS APIs](#). The following diagram shows the software architecture:




Active Robotics Ltd
 1800, 1801, 1802
 1803, 1804, 1805

Active Robotics Ltd
 1800, 1801, 1802, 1803, 1804, 1805



Hardware Specification :

General Assembly [Front]



With Research Robot 2019







Figure 1: 3D Scatter Plot of Data Points

Figure 10: Physical 2D range plot

Figure 11: Physical 2D range plot



Figure 12: Physical 2D range plot



Axis	Range (Degrees)	Range (Meters)
Roll	+90° - 330° - 330°	+0.1407 - 0.1407 - 0.1407
Pitch	+17° - 17°	+0.01 - 0.01 - 0.01
Yaw	+180° - 0° - 180°	+0.0000 - 0.0000 - 0.0000

Open Kinematics: Range of Motion – Total



Joint	Range (Degrees)	Range (Radians)
W1	+179, 0, -179, 359.97	+0, 0.00, -0, 0.00, 0.00
W2	+179.7, 179.7, -179.7, -179.7	+0.00, 0.00, -0.00, -0.00
W3	+179, 0, -179, 359.97	+0, 0.00, -0, 0.00, 0.00
W4	+179.7, 179.7, -179.7, -179.7	+0.00, 0.00, -0.00, -0.00

Joint Position Sensor Based Motion

The resolution for the joint sensors is 14 bits (over 360 degrees) giving $360 / (2^{14}-1) = 0.000023046$ degrees per bit of resolution. All joints have a physical definition of zero, giving a typical accuracy of in the order of ± 0.50 degrees, and not over ± 0.25 degrees accuracy when approaching joint limits. In addition, there may be an absolute zero offset of up to ± 0.50 degrees when the arm is not fully retracted.

Rated Peak Torque

Axis	Peak Torque
W0, W1, W2	100Nm 170N

Positional Accuracy

Type	Accuracy
Wired (No payload) Unloaded (No payload)	$\pm 0.05mm$ 0.1mm

Payload

Type	Accuracy
Max Payload (No End Effector) Max Payload (w/End Effector) (Safety Enabled)	3.5kg 3kg

Single card Configuration

Description	Specification
Processor	Intel Core i7-8070 Processor (8M, 6, 45W)
Memory	16GB DDR4 (RAM)
Hard Drive	256GB Solid State Drive

Camera and Screen Specifications

Description	Specification
Max Resolution	3200 x 2400 pixels
Effective Resolution	2560 x 1920 pixels
Frame Rate	60 frames per second
Focal Length	1.2m
Screen Resolution	3200 x 2400 pixels

Power

Description	Speed Solution
<p>Challenging the application interface</p>	<p>DC-to-3.3VAC converter (Note: the II robot robot has an internal PC, which cannot be powered directly off of 3.3V DC.)</p> <p>Standard 120VAC power is best for sensitive, under-voltage-tolerant devices. Standard power supply cables support 100-240VAC 50/60Hz.</p>
<p>Minimizing energy requirements</p> <p>Electrical Efficiency</p> <p>Power Supply</p>	<p>Optimizing 3.3VAC, 3.3VAC usage per unit is critical.</p> <p>We use standard grade DC and battery power supply for robot of power bus.</p>
<p>Tolerance to surge</p>	<p>We are tolerant to 100% transient voltage surge with respect to power bus.</p>
<p>High speed PCB design</p> <p>Volts range control</p>	<p>High speed design 3.3VDC</p> <p>Surge protection is not needed.</p>